


PRODUCT LIST

Class II Type A2 | CHC-111A2

Airflow



- Room air
- HEPA filtered air
- Contaminated air under negative pressure
- Contaminated air under positive pressure




Model No.		CHC-111A2	
		04	
Dimension (mm)	Width	1340	
	Depth	700	
	Height	1360	
Color		White & Gray	
Door		Tempered safety glass	
Control panel		LCD Display panel	
Sash level		203 mm (8inch)	
External material		SPCC / AL with epoxy powder coated	
Chamber material		STS 304	
Certification		X	
Avg. inflow velocity		0.3 ~ 0.4 m/s	
Avg. Downflow velocity		0.25 ~ 0.35 m/s	
Filter efficiency		HEPA filter > 99.995% at 0.3 μ m	
Noise level (dB)		< 62	
Nominal exhaust volume		312 CMH	
Power supply		220V, 50/60Hz (110V available)	
FL lamp		36W X 1EA	
UV lamp		30W X 1EA	
Power receptacle		U.S type or Euro type (2EA)	
Utility valve		1EA (Vacuum)	
Drain valve		X	
Differential pressure gauge		O	
Duct (Ø)		X	
Hard exhaust connection		Option	
Airtight damper		X	
Exhaust airflow sensor		X	

Support stand (Option)	Support frame (Urethane caster)	Height / mm	
			680
			730

Class II Type A2 | CHC-222A2

Airflow



- Room air
- HEPA filtered air
- Contaminated air under negative pressure
- Contaminated air under positive pressure



Model No.		CHC-222A2			
		03	04	05	06
Dimension (mm)	Width	1060	1385	1690	1995
	Depth	810			
	Height	1590			
Color		White & Gray			
Door		6.38mm Laminated glass			
Control panel		LCD Display panel			
Sash level		203 mm (8inch)			
External material		SPCC / AL with epoxy powder coated			
Chamber material		STS 304			
Certification		-	EN, CE	-	EN, CE
Avg. inflow velocity		0.53 m/s (105fpm)			
Avg. Downflow velocity		0.35 m/s (70fpm)			
Filter efficiency		HEPA filter > 99.995% at 0.3 μ m			
Noise level (dB)		< 64			< 65
Nominal exhaust volume		358 CMH	484 CMH	602 CMH	720 CMH
Power supply		220V, 50/60Hz (110V available)			
FL lamp		36W X 2EA			
UV lamp		30W X 1EA			
Power receptacle		U.S type or Euro type (2EA)			
Utility valve		1EA (Vacuum)			
Drain valve		Option			
Differential pressure gauge		O			
Duct (Ø)		X			
Hard exhaust connection		Option			
Airtight damper		X			
Exhaust airflow sensor		X			

Support stand (Option)	Support frame (Urethane caster)	Height / mm	
			680
			730