



SD-830

Single Drawer Incubator

Open the "Drawer" to new possibilities



Very much like its forerunner, the SD-830 can also be used as a Dry or Humidified incubation unit. Heated from the Top & Bottom of the chamber, the SD-830 promises a stable and well protected environment to preserve/recuperate the temperature inside its chamber.

With a state-of-the-art Gas mixing system placed in the back of the unit, the SD-830 conveniently takes the incoming gas (CO₂- N₂) and provides the necessary mixed levels for your culture.



SPECIFICATIONS

MODEL SD-830

External dimensions (mm)	W310 x D328 x H188	Temp. uniformity	±0.2C
Number of chambers	2	Humidity	90%+
Usable dim. of chamber (mm)	W108 x D180 x H19	Humidity system	Natural evaporation
Weight	10kg	Required gas	CO ₂ /N ₂
Operational environment	24~28C	Power	AC110-115V 2A (Max) 50/60Hz
Temp. accuracy	±0.1C		AC220-240V 1A (Max) 50/60Hz



[Please contact us below]

Email : astec-japan@astec-bio.com TEL : +81-92-935-5666
Or your local ASTEC distributor

www.astec-bio.com/global/